



Bronchoscope Valve Set

Safer, Stronger
and More Sustainable



andorate

The Top 3 Benefits to The Users



Safer with Improved Toggling Mechanism

Working to achieve better patient outcomes by minimising continuous suction that could interrupt the procedure.



Stronger with Advanced Design

Ensuring the suction valve stays intact during installation, procedure and removal.



More Sustainable with Recyclable and Vegetable Ink Printing-based Packaging

Reducing waste with the new bronchoscope valve set by using packaging printed with vegetable-based ink and made with up to 90% recyclable materials.



Scan for more details

Bronchoscope Valves



2-piece Valve kit
(Suction + Biopsy valve)

GAR109

Per box: 50 units
| 315 x 235 x 235 mm | 1.49 kg
Per carton: 200 units
| 490 x 490 x 345 mm | 5.98 kg



Suction Valve

GAR109A

Per box: 50 units
| 315 x 235 x 235 mm | 1.32 kg
Per carton: 200 units
| 490 x 490 x 345 mm | 5.28 kg



Biopsy Valve

GAR021

Per box: 100 units
| 248 x 235 x 174 mm | 0.95 kg
Per carton: 1000 units
| 892 x 482 x 283 mm | 9.60 kg

Compatible with Olympus® 160/180/190/260/290 Series Bronchoscope



We are Global-Asia Health, providing advanced solutions for infection prevention in endoscopy.

GA Health Company Limited
Unit 18, 21/F, Metropole Square
2 On Yiu Street, Shatin, N.T.
Hong Kong
☎ +852 2833 9010
✉ info@GAhealth.com
🌐 GAhealth.com

GA Health Medical Devices Limited
The Glasshouses, Frascati Hall,
Sweetman's Avenue, Blackrock,
Co. Dublin, A94 F9N7, Ireland
✉ info@GAhealth.com
🌐 GAhealth.com