

Endoscope Valve Set

Featuring Champion Design



Innovation for Superior Performance

andorate

Andorate Endoscope Valves

 **Worldwide Patented Design**

 **Advanced Smoothing Technology**

 **Exceptional User Experience**

Multiple Protections

The lubricant-free and pre-sterilised Andorate Endoscope Valves minimise the risk of bacterial interaction with bodily fluids. The single-use design reduces the potential for cross-contamination and eliminates the need for disinfection.

Less Resistance

Designed explicitly for endoscopic applications, the Andorate Endoscope Valves provide an exceptionally smooth user experience by decreasing friction in the valve control mechanisms.

Tactile Feedback

Andorate Endoscope Valves provide precise and instant feedback, enabling the effective and accurate control of air and water during endoscopic procedures.

Reduced Finger Fatigue

Andorate Endoscope Valves require only a light depressing force, significantly reducing discomfort in physician's fingers during procedures.

Eco-friendly Packaging

Packaging of Andorate Endoscope Valve Sets is made from recyclable materials and printed with soy-based ink.



Scan for
product video

Order Information



Endoscope Valve Set

101004

Compatible with Olympus
140/160/180/190/240/260/290/1100/
1200/1500 Series GI Endoscopes

3-piece Set
Air/Water Valve, Suction Valve,
Biopsy Valve

Per box: 50 units
297 x 162 x 187 mm | 1.48 kg
Per carton: 200 units
340 x 312 x 398 mm | 5.94 kg



Endoscope Valve Set with Auxiliary Water Connector

101046

Compatible with Olympus
140/160/180/190/240/260/290/1100/
1200/1500 Series GI Endoscopes

4-piece Set
Air/Water Valve, Suction Valve, Biopsy Valve,
Auxiliary Water Connector

Per box: 50 units
297 x 162 x 187 mm | 1.77 kg
Per carton: 200 units
340 x 312 x 398 mm | 7.10 kg



Endoscope Valve Set with Air/Water Cleaning Adapter

101047

Compatible with Olympus
140/160/180/190/240/260/290/1100/
1200/1500 Series GI Endoscopes

4-piece Set
Air/Water Valve, Suction Valve, Biopsy Valve,
Air/Water Cleaning Adapter

Per box: 20 units
285 x 175 x 220 mm | 0.86 kg
Per carton: 80 units
300 x 370 x 470 mm | 4.10 kg



Endoscope Valve Set with Auxiliary Water Connector & Air/Water Cleaning Adapter

101048

Compatible with Olympus
140/160/180/190/240/260/290/1100/
1200/1500 Series GI Endoscopes

5-piece Set
Air/Water Valve, Suction Valve, Biopsy Valve,
Auxiliary Water Connector,
Air/Water Cleaning Adapter

Per box: 20 units
285 x 175 x 220 mm | 0.90 kg
Per carton: 80 units
300 x 370 x 470 mm | 4.27 kg



Endoscope Combo Valve Set

101050

Compatible with Olympus
140/160/180/190/240/260/290/1100/
1200/1500 Series GI Endoscopes

6-piece Set
Air/Water Valve, Suction Valve, Biopsy Valve,
Auxiliary Water Connector, Air/Water Cleaning
Adapter, Biopsy Port Cleaning Adapter

Per box: 20 units
285 x 175 x 220 mm | 1.37 kg
Per carton: 80 units
300 x 370 x 470 mm | 6.14 kg

● Compatible with Olympus endoscopes



For single-patient use

Order Information



Endoscope Valve Set

101087

Compatible with Fujifilm
700 series GI Endoscopes

3-piece Set
Air/Water Valve, Suction Valve,
Biopsy Valve

Per box: 50 units
297 x 162 x 187 mm | 1.66 kg
Per carton: 200 units
340 x 312 x 398 mm | 6.66 kg



Endoscope Valve Set with Auxiliary Water Connector

101092

Compatible with Fujifilm
700 series GI Endoscopes

4-piece Set
Air/Water Valve, Suction Valve,
Biopsy Valve, Auxiliary Water
Connector

Per box: 50 units
297 x 162 x 187 mm | 1.73 kg
Per carton: 200 units
340 x 312 x 398 mm | 6.92 kg



Endoscope Valve Set with Air/Water Cleaning Adapter

101093

Compatible with Fujifilm
700 series GI Endoscopes

4-piece Set
Air/Water Valve, Suction Valve,
Biopsy Valve, Air/Water
Cleaning Adapter

Per box: 20 units
285 x 175 x 220 mm | 0.83 kg
Per carton: 80 units
300 x 370 x 470 mm | 4.00 kg



Endoscope Valve Set with Auxiliary Water Connector & Air/Water Cleaning Adapter

101094

Compatible with Fujifilm
700 series GI Endoscopes

5-piece Set
Air/Water Valve, Suction Valve,
Biopsy Valve, Auxiliary Water Connector,
Air/Water Cleaning Adapter

Per box: 20 units
285 x 175 x 220 mm | 0.87 kg
Per carton: 80 units
300 x 370 x 470 mm | 4.14 kg



Endoscope Combo Valve Set

101096

Compatible with Fujifilm
700 series GI Endoscopes

6-piece Set
Air/Water Valve, Suction Valve,
Biopsy Valve, Auxiliary Water Connector,
Air/Water Cleaning Adapter,
Biopsy Port Cleaning Adapter

Per box: 20 units
285 x 175 x 220 mm | 1.33 kg
Per carton: 80 units
300 x 370 x 470 mm | 5.99 kg

● Compatible with Fujifilm endoscopes



For single-patient use

Order Information



Endoscope Valve Set

101037

Compatible with Pentax
90/i10/i10c Series GI Endoscopes

3-piece Set
Air/Water Valve, Suction Valve,
Biopsy Valve

Per box: 50 units
297 x 162 x 187 mm | 1.65 kg
Per carton: 200 units
340 x 312 x 398 mm | 6.62 kg



Endoscope Valve Set with Auxiliary Water Connector

101081

Compatible with Pentax
90/i10/i10c Series GI Endoscopes

4-piece Set
Air/Water Valve, Suction Valve,
Biopsy Valve, Auxiliary Water
Connector

Per box: 50 units
297 x 162 x 187 mm | 1.73 kg
Per carton: 200 units
340 x 312 x 398 mm | 6.92 kg

• Compatible with Pentax endoscopes



For single-patient use

Please contact Andorate authorised distributor in your country or email us at sales@GAhealth.com for ordering details.



Scan and Contact Us



andorate



We are Global-Asia Health,
providing advanced solutions
for infection prevention in
endoscopy.

GA Health Company Limited

Unit 1511-13 & 15, 15/F,
Tower 2, Grand Central Plaza,
138 Shatin Rural Committee Road,
Shatin, N.T., Hong Kong

 +852 2833 9010

 info@GAhealth.com

 GAhealth.com

GA Health Medical Devices Limited

The Glasshouses, Frascati Hall,
Sweetman's Avenue, Blackrock,
Co. Dublin, A94 F9N7, Ireland

 info@GAhealth.com

 GAhealth.com