

Andorate **CHAMPION** Endoscope Valves

Innovative Technology for Superior Performance



Less Resistance • Tactile Feedback • Reduced Finger Fatigue

andorate

Andorate **CHAMPION** Endoscope Valves

 **Worldwide Patented Design**

 **Advanced Smoothing Technology**

 **Exceptional User Experience**



Multiple Protections

The lubricant-free and pre-sterilised CHAMPION Endoscope Valve minimises the risk of bacterial interaction with bodily fluids. The single-use design reduces the potential for cross-contamination and eliminates the need for disinfection.



Less Resistance

Designed explicitly for endoscopic applications, the CHAMPION Endoscope Valve provides an exceptionally smooth user experience by decreasing friction in the valve control mechanisms.



Tactile Feedback

CHAMPION Endoscope Valve provides precise and instant feedback, enabling the effective and accurate control of air and water during endoscopic procedures.



Reduced Finger Fatigue

CHAMPION Endoscope Valve requires only a light depressing force, significantly reducing discomfort in physician's fingers during procedures.

Eco-friendly Packaging

Packaging of CHAMPION is made from recyclable materials and printed with soy-based ink.



Scan for
CHAMPION video

Order Information



CHAMPION Valve Set

101004

Compatible with Olympus
140/160/180/190/240/260/290/1100/
1200/1500 Series GI Endoscopes

3-piece Set
Air-Water Valve, Suction Valve,
Biopsy Valve

Per box: 50 units
297 x 162 x 187 mm | 1.48 kg
Per carton: 200 units
340 x 312 x 398 mm | 5.94 kg



CHAMPION Valve Set with Auxiliary Water Connector

101046

Compatible with Olympus
140/160/180/190/240/260/290/1100/
1200/1500 Series GI Endoscopes

4-piece Set
Air-Water Valve, Suction Valve, Biopsy Valve,
Auxiliary Water Connector

Per box: 50 units
297 x 162 x 187 mm | 1.77 kg
Per carton: 200 units
340 x 312 x 398 mm | 7.10 kg



CHAMPION Valve Set with Air-Water Cleaning Adapter

101047

Compatible with Olympus
140/160/180/190/240/260/290/1100/
1200/1500 Series GI Endoscopes

4-piece Set
Air-Water Valve, Suction Valve, Biopsy Valve,
Air-Water Cleaning Adapter

Per box: 20 units
285 x 175 x 220 mm | 0.86 kg
Per carton: 80 units
300 x 370 x 470 mm | 4.10 kg



CHAMPION Valve Set with Auxiliary Water Connector & Air-Water Cleaning Adapter

101048

Compatible with Olympus
140/160/180/190/240/260/290/1100/
1200/1500 Series GI Endoscopes

5-piece Set
Air-Water Valve, Suction Valve, Biopsy Valve,
Auxiliary Water Connector,
Air-Water Cleaning Adapter

Per box: 20 units
285 x 175 x 220 mm | 0.90 kg
Per carton: 80 units
300 x 370 x 470 mm | 4.27 kg



CHAMPION Combo Valve Set

101050

Compatible with Olympus
140/160/180/190/240/260/290/1100/
1200/1500 Series GI Endoscopes

6-piece Set
Air-Water Valve, Suction Valve, Biopsy Valve,
Auxiliary Water Connector, Air-Water Cleaning
Adapter, Biopsy Port Cleaning Adapter

Per box: 20 units
285 x 175 x 220 mm | 1.37 kg
Per carton: 80 units
300 x 370 x 470 mm | 6.14 kg

● Compatible with Olympus endoscopes



For single-patient use

Order Information



CHAMPION Valve Set

101087

Compatible with Fujifilm
700 series GI Endoscopes

3-piece Set
Air-Water Valve, Suction Valve,
Biopsy Valve

Per box: 50 units
297 x 162 x 187 mm | 1.66 kg
Per carton: 200 units
340 x 312 x 398 mm | 6.66 kg



CHAMPION Valve Set with Auxiliary Water Connector

101092

Compatible with Fujifilm
700 series GI Endoscopes

4-piece Set
Air-Water Valve, Suction Valve,
Biopsy Valve, Auxiliary Water
Connector

Per box: 50 units
297 x 162 x 187 mm | 1.73 kg
Per carton: 200 units
340 x 312 x 398 mm | 6.92 kg



CHAMPION Valve Set with Air-Water Cleaning Adapter

101093

Compatible with Fujifilm
700 series GI Endoscopes

4-piece Set
Air-Water Valve, Suction Valve,
Biopsy Valve, Air-Water
Cleaning Adapter

Per box: 20 units
285 x 175 x 220 mm | 0.83 kg
Per carton: 80 units
300 x 370 x 470 mm | 4.00 kg



CHAMPION Valve Set with Auxiliary Water Connector & Air-Water Cleaning Adapter

101094

Compatible with Fujifilm
700 series GI Endoscopes

5-piece Set
Air-Water Valve, Suction Valve,
Biopsy Valve, Auxiliary Water Connector,
Air-water Cleaning Adapter

Per box: 20 units
285 x 175 x 220 mm | 0.87 kg
Per carton: 80 units
300 x 370 x 470 mm | 4.14 kg



CHAMPION Combo Valve Set

101096

Compatible with Fujifilm
700 series GI Endoscopes

6-piece Set
Air-Water Valve, Suction Valve,
Biopsy Valve, Auxiliary Water Connector,
Air-Water Cleaning Adapter,
Biopsy Port Cleaning Adapter

Per box: 20 units
285 x 175 x 220 mm | 1.33 kg
Per carton: 80 units
300 x 370 x 470 mm | 5.99 kg

● Compatible with Fujifilm endoscopes



For single-patient use

Order Information



CHAMPION Valve Set

101037

Compatible with Pentax
90/i10/i10c Series GI Endoscopes

3-piece Set
Air-Water Valve, Suction Valve,
Biopsy Valve

Per box: 50 units
297 x 162 x 187 mm | 1.65 kg
Per carton: 200 units
340 x 312 x 398 mm | 6.62 kg



CHAMPION Valve Set with Auxiliary Water Connector

101081

Compatible with Pentax
90/i10/i10c Series GI Endoscopes

4-piece Set
Air-Water Valve, Suction Valve,
Biopsy Valve, Auxiliary Water
Connector

Per box: 50 units
297 x 162 x 187 mm | 1.73 kg
Per carton: 200 units
340 x 312 x 398 mm | 6.92 kg

● Compatible with Pentax endoscopes



For single-patient use

Please contact Andorate authorised distributor in your country or email us at sales@GAhealth.com for ordering details.



Scan and Contact Us



andorate



We are Global-Asia Health,
providing advanced solutions
for infection prevention in
endoscopy.

GA Health Company Limited

Unit 1511-13 & 15, 15/F,
Tower 2, Grand Central Plaza,
138 Shatin Rural Committee Road,
Shatin, N.T., Hong Kong

☎ +852 2833 9010

✉ info@GAhealth.com

🌐 GAhealth.com

GA Health Medical Devices Limited

The Glasshouses Blackrock, Frascati Hall,
Sweetman's Ave, Newtown Blackrock,
Blackrock, Co. Dublin, A94 F9N7, Ireland

✉ info@GAhealth.com

🌐 GAhealth.com