

Andorate **PROTEGO** Air-Water System

Preventing the Contaminated Fluid Travelling
Back to the Water Bottle

Defense • Safety • Simplicity



andorate

Andorate **PROTEGO Air-Water System**

How to Reduce the Risk of Backflow During a Gastroenterology Procedure?



DEFENSE

Worldwide Patented Backflow Prevention Technology

PROTEGO Air-Water System functions as a shield to prevent contaminated fluid from travelling back to the water bottle, improving patient and user safety.



SAFETY

Efficient Care Solution


PROTEGO eliminates the need to replace the tubing and water bottles for high-risk patients, saving money and reducing waste.



SIMPLICITY

Designed for User Convenience

PROTEGO Air-Water System features colour-coded tubing and connector for easy identification and a one-click open and close mechanism for straightforward operation.



Scan for
PROTEGO video



PROTEGO Air-Water Tubing
for 24-hour multi-patient use

PROTEGO Air-Water Connector
for single-patient use

Order Information



24h

PROTEGO Air-Water Bottle Tubing

202001

Compatible with Olympus
160/180/190/240/260/290/1100/1200/
1500 Series GI Endoscopes

For use with the
PROTEGO Air-Water Connector

Per box: 20 units
380 x 300 x 150 mm | 2.85 kg
Per carton: 40 units
425 x 350 x 310 mm | 5.70 kg



24h

PROTEGO Air-Water Bottle Tubing with CO₂ Extend

202001CO2EXT

Compatible with Olympus
160/180/190/240/260/290/1100/1200/
1500 Series GI Endoscopes

For use with the
PROTEGO Air-Water Connector

Per box: 10 units
297 x 162 x 202 mm | 2.10 kg
Per carton: 40 units
425 x 350 x 310 mm | 8.50 kg



2

PROTEGO Air-Water Connector

201001

Compatible with Olympus
160/180/190/240/260/290/1100/1200/
1500 Series GI Endoscopes

For use with the
PROTEGO Air-Water Bottle Tubing

Per box: 100 units
380 x 300 x 150 mm | 1.40 kg
Per carton: 200 units
425 x 350 x 310 mm | 3.00 kg



24h

PROTEGO Air-Water Bottle Tubing

202002

Compatible with Fujifilm
500/600/700 series GI Endoscopes

For use with the
PROTEGO Air-Water Connector

Per box: 20 units
380 x 300 x 150 mm | 2.85 kg
Per carton: 40 units
425 x 350 x 310 mm | 5.70 kg



24h

PROTEGO Air-Water Bottle Tubing with CO₂ Extend

202002CO2EXT

Compatible with Fujifilm
500/600/700 series GI Endoscopes

For use with the
PROTEGO Air-Water Connector

Per box: 10 units
297 x 162 x 202 mm | 1.40 kg
Per carton: 40 units
425 x 350 x 310 mm | 6.00 kg



2

PROTEGO Air-Water Connector

201002

Compatible with Fujifilm
700 series GI Endoscopes

For use with the
PROTEGO Air-Water Bottle Tubing

Per box: 100 units
315 x 235 x 235 mm | 2.00 kg
Per carton: 400 units
490 x 490 x 345 mm | 8.43 kg



2

PROTEGO Air-Water Connector

201003

Compatible with Fujifilm
500/600 series GI Endoscopes

For use with the
PROTEGO Air-Water Bottle Tubing

Per box: 100 units
315 x 235 x 235 mm | 2.00 kg
Per carton: 400 units
490 x 490 x 345 mm | 8.43 kg

● Compatible with Olympus endoscopes ● Compatible with Fujifilm endoscopes



For single-patient use



For 24-hour use

Please contact Andorate authorised distributor in your country or email us at sales@GAhealth.com for ordering details.



Scan and Contact Us

andorate



We are Global-Asia Health,
providing advanced solutions
for infection prevention in
endoscopy.

GA Health Company Limited

Unit 1511-13 & 15, 15/F,
Tower 2, Grand Central Plaza,
138 Shatin Rural Committee Road,
Shatin, N.T., Hong Kong

☎ +852 2833 9010

✉ info@GAhealth.com

🌐 GAhealth.com

GA Health Medical Devices Limited

The Glasshouses, Frascati Hall,
Sweetman's Avenue, Blackrock,
Co. Dublin, A94 F9N7, Ireland

✉ info@GAhealth.com

🌐 GAhealth.com